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Thank you for purchasing this dental product.

Carefully read through these instructions before each use and store them somewhere that is easy for the user and for the relevant qualified personnel to access.

Carefully read the warning information indicated by this warning symbol. Improper use of the products can lead to serious injury to the patient, user or any third parties. To ensure that the condition and functionality is faultless for the intended use, the specifications of the following instructions must be complied with. Please note that handling the implements improperly can negatively affect their service life and safety.

The medical devices are not delivered in a sterile state and must be prepared and sterilised by the user according to the following instructions before the first and before any subsequent use.

## 1 Product details:

These user instructions are applicable for the products specified below:

|     |                                                                                                                                                                                                                                                                                                                                                                                                                                                                             |                     |
|-----|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|---------------------|
| REF | 301-01-R, 301-04-R, 301-11-R, 303-02-R, 303-10-R, 303-16-R, 304-02-R, 304-11-R, 304-26-R, 308-01-R, 308-04-R, 308-09-R, 313-01-R, 313-05-R, 320-02-R, 320-04-R, 320-08-R.                                                                                                                                                                                                                                                                                                   | (R = Stainless)     |
|     | 301-01-N, 301-03-N, 301-04-N, 301-06-N, 301-07-N, 301-12-N, 303-01-N, 303-03-N, 303-05-N, 303-06-N, 303-09-N, 303-15-N, 303-17-N, 303-18-N, 304-01-N, 304-03-N, 304-04-N, 304-05-N, 304-06-N, 304-12-N, 304-19-N, 304-21-N, 304-24-N, 304-30-N, 304-31-N, 308-01-N, 308-02-N, 308-04-N, 308-05-N, 308-07-N, 308-08-N, 308-09-N, 310-01-N, 310-03-N, 310-06-N, 310-09-N, 313-01-N, 313-05-N, 320-01-N, 320-03-N, 320-05-N, 320-06-N, 320-07-N, 320-08-N, 320-09-N, 320-11-N. | (N = Nickel-plated) |

## 2 Purpose:

**MD** This product is a medical device.

Mandrels are used for stretching elastic polishers, paper discs, nylon discs, diamond-coated disks and brushes for application in the field of dentistry. The products have attachments conforming to DIN EN ISO 1797 – Type 2 (HP shank Ø 2.35 mm) or Type 4 (HP shank Ø 3.00 mm).

### 2.1 Indication

Corrosion-resistant, stainless steel (R) mandrels are intended for use in invasive and non-invasive applications.

Nickel-plated (N) mandrels are intended for use in non-invasive applications.

### 2.2 Contraindication

This product must not be used in case of an intolerance or allergy to standard surgical stainless steel.

### 2.3 Intended use

The implements should only be used for their intended purpose in the field of dentistry by trained and qualified personnel. The treating doctor, buyer or user is responsible for selecting the appropriate implements for certain applications or for operative use. Josef Ganter is not liable for any direct damage or consequential damage caused by improper use, handling, preparation or maintenance.

### 2.4 Warnings → Use / Inspection / Usability / Safety

- May only be used by an adequately qualified person.
- Always wear protective glasses.
- Stretched tools should not have any unbalance since the mandrel could break. Tools which do not run smoothly must be adjusted at a low rotational speed before use.
- Never exceed the recommended speed for the installed article (grinding tool, polisher, etc.).
- Do not continue to use damaged mandrels (e.g.: damaged plating, screws can no longer be sufficiently tightened, mandrel does no longer run true). Perform a visual inspection before each use.
- Mandrels with screws are only intended for use in clockwise rotation (exception: mandrels with left-hand thread – only for left-hand rotation).
- Only use mandrels with screws with washers. Avoid excessively tightening the screws (risk of breaking).
- Avoid sudden changes in speed.
- Work with a moderate contact pressure (max. 3 N).
- HP mandrels must always be clamped along the entire length of the shank.
- Larger or heavier tools, such as those mentioned under point “7 Combination with other products and permissible speeds”, should not be attached to the mandrels.

## 3 Preparation

### 3.1 Warnings → preparation

- The preparation specification is only valid for the mandrel without combination products such as grinding tools, polishers, etc.. Combination products must be removed for preparation and prepared according to the manufacturer's instructions.
- The preparation may only be performed by adequately qualified persons.
- The water used must be at least of drinking water quality. (see specifications in the individual preparation steps)
- In these preparation instructions, the cleaning and disinfecting agents used for the efficiency verification will be specified. If you use an alternative cleaning and / or disinfecting agent, it must be listed by the RKI [Robert Koch Institute] or the VAH [German Association for Applied Hygiene] and must be compatible with the materials. The pH value must be between 4.5 and 10.
- The preparator is responsible for achieving the desired results when actually performing the preparation in the preparation facility with the equipment, material and persons used. Therefore, as a general rule, the validation and routine monitoring of the procedure and the equipment used is required.
- A mechanical cleaning and disinfection process (cleaning and disinfection unit) is currently not permissible.

## 3.2 Transport / site of use – preparation

The first step in preparing a product correctly starts immediately after it has been used on a patient.

Heavy contamination, residues of fillings, disinfection agents and other medicinal products should be removed before the implements are stored away.

- **Dry removal (humidified, closed system) is to be preferred whenever and wherever possible.**
- **As a general rule, surface drying of certain residues which are left after use is to be avoided!**
- **Long waiting periods before the preparation, e.g. overnight or over the weekend, are to be avoided with both types of removal (<6 hours).**

## 3.3 Preparation for decontamination

The products are to be, insofar as possible, disassembled before the following preparation steps and the further preparation steps are to be carried out when they are open. Any traces of cleaning agents are to be avoided. The products must be prepared in suitable sieve baskets or drainage bowls (choose the size depending on the product). A minimum distance should be maintained between the products in the cleaning basket. Overlapping is to be avoided to prevent any damage to the product during the cleaning process.

## 3.4 Pre-cleaning

Pre-cleaning must always be performed regardless of the following cleaning method. Rinse the products under cold municipal water (drinking water quality, <40°C) until all visible contamination has been removed. Any dirt still adhering to the product must be removed with a soft brush. Hollow spaces and lumens must be intensively (>30 seconds) rinsed out using a water pistol (or similar) with cold municipal water (drinking water quality <40°C). Moving parts must be moved whilst they are being cleaned.

## 3.5 Cleaning / disinfection

### 3.5.1 Manual cleaning and disinfection process

1. Place products in a cleaning and disinfection agent (2% *Stabimed*) in an ultrasound bath for approx. 15 minutes and at a frequency of 35 kHz, at room temperature.
2. Thoroughly clean the product with a soft brush afterwards. If there are any hollow spaces and lumens, intensively (>30 seconds) rinse them out using a water pistol (or similar).
3. Rinse the products (complete rinsing of the inside, outside and hollow spaces) in demineralised water for >60 seconds.

### 3.5.2 Manual drying

Dry manually with a lint-free, single-use cloth. To avoid leaving any water in hollow spaces, blow them out with sterile, oil-free pressurised air. The product should not be heated to over 140°C.

## 3.6 Sterilisation – steam sterilisation – fractionated vacuum procedure

*Evidence of the fundamental suitability of the implements for an effective sterilisation was provided by an independent and accredited testing laboratory that is recognised by the ZLG [Central Authority of the Länder for Health Protection with regard to Medicinal Products and Medical Devices] (§15 (5) MGP [Medical Products Directive]) under use of an EHS3870 pre- and post-vacuum autoclave (Tuttnauer Europe B.V., Breda) and the sterilised packaging used. For this, the procedure described above was used. These specifications must be complied with.*

*2 vacuum cycles | 134°C / 273°F | ≥ 2.5 minutes stop time | Dry in a vacuum for at least 15 minutes<sup>3</sup>*

### 3.6.1 Packaging

The sterilisation of the products must be done in suitable sterilisation packaging. The manufacturer's verification was done using doubled sterilisation packaging (hospital standard); that means that the single suitable sterilised foil packaging can also be used.

**Flash sterilisation and the sterilisation of unpackaged implements is absolutely prohibited!**

### 3.6.2 Sterilisation process – Fractionated vacuum procedure

Only the specified sterilisation procedures may be used for sterilisation.

Other sterilisation procedures are not permitted and their efficiency must be certified by the user / preparator themselves.

- **Fractionated vacuum procedure<sup>1,2</sup> (with sufficient product drying<sup>3</sup>)**
- Steam steriliser conforming to DIN EN 13060/DIN EN 285 and ANSI AAMI ST 79 (FDA clearance in the USA)
- Validated in accordance with DIN EN ISO 17665 (valid IQ/OQ (commissioning) and product-specific performance evaluation (PQ))
- Maximum sterilisation temperature 134°C (273°F) including the tolerance margin in accordance with DIN EN ISO 17665
- Sterilisation time (time exposed to the sterilisation temperature)

## 4 Storage

After the sterilisation, the products must be stored dust-free and dry in the sterilisation packaging.

## 5 Product life

The service life of the disc mandrel ends after 50 fastening cycles (number of assemblies or disassemblies of the grinding tool).

## 6 Repairs

Josef Ganter GmbH provides a repair service for your implements (insofar as this is possible). If the implements are repaired by companies or persons which are not expressly authorised to do so by Josef Ganter, all guarantees shall become null and void.

<sup>1</sup> At least 3 vacuum stages

<sup>2</sup> The less effective gravitation procedure may only be used if the fractionated vacuum procedure is not available as an option, as it requires significantly longer sterilisation periods which must be determined and validated under the responsibility of the user, and specifically for their implements, devices, procedures and parameters.

<sup>3</sup> The actual required product drying time depends directly on the parameters, which are solely the responsibility of the user (loading configuration and density, sterilisation status,...) and must therefore be determined by them. As a general rule, drying times should not be less than 15 minutes.

## 7 Combination with other products and permissible speeds

| Maximum permissible speed in revolutions per minute                                                        |                                        |                                                           |                                  |                                  |                                  |                                  |                               |                                 |                    |                                                  |
|------------------------------------------------------------------------------------------------------------|----------------------------------------|-----------------------------------------------------------|----------------------------------|----------------------------------|----------------------------------|----------------------------------|-------------------------------|---------------------------------|--------------------|--------------------------------------------------|
| Article Ref. no.                                                                                           | Maximum programming length of the tool | Shank connection for straight and contra-angle handpieces | Discs with max. Ø 10 x 2 (< 1 g) | Discs with max. Ø 22 x 1 (< 1 g) | Discs with max. Ø 22 x 3 (< 3 g) | Discs with max. Ø 38 x 1 (< 3 g) | Moore-Discs max. Ø 15 (< 1 g) | Polishing pins Length max. 23mm | Sand-paper stripes | Felts and polisher with max. Ø7x 25 long (< 3 g) |
| 303-01<br>303-02<br>303-03<br>303-05<br>304-01<br>304-02<br>304-04<br>304-05<br>304-06<br>304-26<br>304-30 | 25mm                                   | HP   2,35 mm                                              | 20.000                           | 20.000                           | 10.000                           | -                                | -                             | -                               | -                  | -                                                |
| 303-06<br>304-03<br>304-19                                                                                 | 25mm                                   | HP   3,00 mm                                              | 20.000                           | 20.000                           | 10.000                           | -                                | -                             | -                               | -                  | -                                                |
| 303-09<br>303-10                                                                                           | 25mm                                   | HP   2,35 mm                                              | 20.000                           | -                                | -                                | -                                | -                             | -                               | -                  | -                                                |
| 303-15<br>303-16<br>303-17<br>304-11<br>304-12<br>304-24                                                   | 25mm                                   | HP   2,35 mm                                              | 15.000                           | 15.000                           | 10.000                           | 5.000                            | -                             | -                               | -                  | -                                                |
| 303-18<br>304-21<br>304-31                                                                                 | 25mm                                   | HP   3,00 mm                                              | 15.000                           | 15.000                           | 10.000                           | 5.000                            | -                             | -                               | -                  | -                                                |
| 313-01                                                                                                     | 25mm                                   | HP   2,35 mm                                              | -                                | -                                | -                                | -                                | 20.000                        | -                               | -                  | -                                                |
| 313-05                                                                                                     | 25mm                                   | HP   2,35 mm                                              | -                                | -                                | -                                | -                                | 10.000                        | -                               | -                  | -                                                |
| 308-01<br>308-04<br>308-07<br>308-08<br>308-09                                                             | 30mm                                   | HP   2,35 mm                                              | -                                | -                                | -                                | -                                | -                             | -                               | 15.000             | -                                                |
| 308-02<br>308-05                                                                                           | 30mm                                   | HP   3,00 mm                                              | -                                | -                                | -                                | -                                | -                             | -                               | 15.000             | -                                                |
| 320-01<br>320-02<br>320-03<br>320-04                                                                       | 25mm                                   | HP   2,35 mm                                              | -                                | -                                | -                                | -                                | -                             | 20.000                          | -                  | -                                                |
| 320-05                                                                                                     | 25mm                                   | HP   2,35 mm                                              | -                                | -                                | -                                | -                                | -                             | 15.000                          | -                  | -                                                |
| 320-06<br>320-07<br>320-08<br>320-09<br>320-11                                                             | 25mm                                   | HP   2,35 mm                                              | -                                | -                                | -                                | -                                | -                             | 10.000                          | -                  | -                                                |
| 301-01<br>301-04<br>301-07<br>301-11<br>301-12                                                             | 25mm                                   | HP   2,35 mm                                              | -                                | -                                | -                                | -                                | -                             | -                               | -                  | 10.000                                           |
| 301-03<br>301-06                                                                                           | 25mm                                   | HP   3,00 mm                                              | -                                | -                                | -                                | -                                | -                             | -                               | -                  | 10.000                                           |
| 310-01                                                                                                     | 25mm                                   | HP   2,35 mm                                              | -                                | -                                | -                                | -                                | -                             | -                               | -                  | 10.000                                           |
| 310-03<br>313-06<br>313-09                                                                                 | 25mm                                   | HP   3,00 mm                                              | -                                | -                                | -                                | -                                | -                             | -                               | -                  | 10.000                                           |

| Symbols used                                                                        |                                                        |
|-------------------------------------------------------------------------------------|--------------------------------------------------------|
| Symbol                                                                              | Description / Symbol for ...                           |
|  | ... Manufacturer                                       |
|  | ... Article number                                     |
|  | ... Production lot number – batch                      |
|  | ... Observe user instructions                          |
|  | ... CE conformity with the number of the notified body |
|  | ... Is a medical device                                |
|  | ... Non-sterile                                        |
|  | ... Warning                                            |

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**Josef Ganter Feinmechanik GmbH**