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info@josefganter.de

www.josefganter.de

Thank you for purchasing this dental product.

Carefully read through these instructions before each use and store them somewhere that is easy for the user and for the relevant qualified personnel to access.

Carefully read the warning information indicated by this warning symbol. Improper use of the products can lead to serious injury to the patient, user or any third parties. To ensure that the condition and functionality is faultless for the intended use, the specifications of the following instructions must be complied with. Please note that handling the implements improperly can negatively affect their service life and safety.

The medical devices are not delivered in a sterile state and must be prepared and sterilised by the user according to the following instructions before the first and before any subsequent use.

1 Product details

These user instructions are applicable for the products specified below:

320-13-R, 320-16-R, 320-17-R

(Stainless **(R)** – corrosion-resistant stainless steel)

2 Purpose

MD This product is a medical device.

Drill extensions are for extending articles for the field of dentistry. The products have attachments conforming to DIN EN ISO 1797 – Type 1 (RA-shank 2.35 mm).

2.1 Indication and combination

Drill extensions are intended for use as an adapter for the extension of tools with a shank conforming to DIN EN ISO 1797 (Type 1 – RA-shank) and the standardised connector for the contra-angle handpiece, to allow for better accessibility to hard-to-reach areas during dental procedures.

2.2 Contraindication

This product must not be used in case of an intolerance or allergy to standard surgical stainless steel.

2.3 Intended use

The implements should only be used for their intended purpose in the field of dentistry by trained and qualified personnel. The treating doctor, buyer or user is responsible for selecting the appropriate implements for certain applications or for operative use. Josef Ganter is not liable for any direct damage or consequential damage caused by improper use, handling, preparation or maintenance.

2.4 Warnings → Use / Inspection / Usability / Safety

- May only be used by an adequately qualified person.
- Always wear protective glasses.
- Stretched tools should not have any unbalance since the mandrel could break. Tools which do not run smoothly must be adjusted at a low rotational speed before use.
- Never exceed the recommended speed for the installed article (grinding tool, polisher, etc.).
- Damaged extensions:
e.g.: locking mechanism on the tool or contra-angle handpiece no longer functions, extension does no longer run true → do not continue to use. Perform a visual inspection and functional check before each use.
- Only use intact and unspoiled implements / tools with an RA connection.
- Avoid sudden changes in speed.
- Work with a moderate contact pressure (max. 3 N).
- Do not exceed a speed of 10,000 revolutions per minute.
- RA shanks must always be clamped along the entire length of the shank.

3 Preparation

3.1 Warnings → preparation

- The preparation may only be performed by adequately qualified persons.
- The water used must be at least of drinking water quality. (see specifications in the individual preparation steps).
- In these preparation instructions, the cleaning and disinfecting agents used for the efficiency verification will be specified. If you use an alternative cleaning and / or disinfecting agent, it must be listed by the *RKI* [Robert Koch Institute] or the *VAH* [German Association for Applied Hygiene] and must be compatible with the materials. The pH value must be between 4.5 and 10.
- The preparator is responsible for achieving the desired results when actually performing the preparation in the preparation facility with the equipment, material and persons used. Therefore, as a general rule, the validation and routine monitoring of the procedure and the equipment used is required.

3.2 Transport / site of use – preparation

The first step in preparing a product correctly starts immediately after it has been used on a patient.

Heavy contamination, residues of fillings, disinfection agents and other medicinal products should be removed before the implements are stored away.

- Dry removal (humidified, closed system) is to be preferred whenever and wherever possible.
- As a general rule, surface drying of certain residues which are left after use is to be avoided!
- Long waiting periods before the preparation, e.g. overnight or over the weekend, are to be avoided with both types of removal (<6 hours).

3.3 Preparation for decontamination

The products are to be separated from all tools and contra-angle handpieces before continuing with the following preparation steps. Any traces of cleaning agents are to be avoided. The products must be prepared in suitable sieve baskets or drainage bowls (choose the size depending on the product). A minimum distance should be maintained between the products in the cleaning basket. Overlapping is to be avoided to prevent any damage to the product during the cleaning process.

3.4 Pre-cleaning

Pre-cleaning must always be performed regardless of the following cleaning method. Rinse the products under cold municipal water (drinking water quality, <40°C) until all visible contamination has been removed. Any dirt still adhering to the product must be removed with a soft brush. Hollow spaces and lumens must be intensively (>30 seconds) rinsed out using a water pistol (or similar) with cold municipal water (drinking water quality <40°C).

3.5 Cleaning / disinfection

As a general rule mechanical cleaning is preferred to manual cleaning!

3.5.1 Manual process**3.5.1.1** **Manual cleaning**

1. Place products in an alkaline cleaning agent (0.5% Neodisher® MediClean) in an ultrasonic bath for approx. 10 minutes. Do not exceed the maximum temperature of 40°C. Here, the instructions provided by the cleaning agent manufacturer must be followed.
2. Thoroughly clean the product with a soft brush afterwards. If there are any hollow spaces and lumens, intensively (>30 seconds) rinse them out using a water pistol (or similar).
3. Rinse the product under running municipal water (drinking water quality) to remove the cleaning agent (>15 seconds).

3.5.1.2 **Manual disinfection**

1. Immerse the product in an RKI- or VAH-listed disinfecting agent. Here, the instructions provided by the disinfecting agent manufacturer must be followed. It must be ensured that the disinfecting agent really reaches all areas of the product (move the parts around in the disinfection bath and, if necessary, rinse hidden areas using a syringe – without a cannula – with the disinfecting agent).
2. The efficiency verification for the process was done using the disinfecting agent: 3% Korsolex plus (Bode Chemie, Hamburg) 15 minutes.
3. Rinse the products (complete rinsing of the inside, outside and hollow spaces) in demineralised water for >60 seconds.

3.5.1.3 **Manual drying**

1. Dry manually with a lint-free, single-use cloth. To avoid leaving any water in hollow spaces, blow them out with sterile, oil-free pressurised air. The product should not be heated to over 140°C.

3.5.2 Mechanical process – thermal disinfection

Evidence of the fundamental suitability of the implements for an effective mechanical cleaning and disinfection was provided by an independent and accredited testing laboratory that is recognised by the ZLG [Central Authority of the Länder for Health Protection with regard to Medicinal Products and Medical Devices] (§15 (5) MGP [Medical Products Directive]) under use of the Miele G7835 CD cleaning and disinfection unit (thermal disinfection, Miele & Cie. KG, Gütersloh) and the pre-cleaning and cleaning agents neodisher® mediClean (Dr. Weigert GmbH & Co. KG, Hamburg). For this, the procedure described above was used. These specifications must be complied with.

3.5.2.1 **Cleaning and disinfection unit and media**

When choosing a cleaning and disinfection unit, it should be ensured:

- that the effectiveness of the unit has been verified (e.g. DGHM [German Association for Hygiene and Microbiology] or FDA [Food and Drug Administration] approval / clearance / registration or CE label in accordance with DIN EN ISO 15883),
- that where possible a tested thermal disinfection programme (A₀ value > 3000 or – with older devices – at least 5 minutes at 90°C / 194°F) is used (with chemical disinfection there is a risk of disinfecting agent residue on the implements),
- that the programme used is suitable for the implements and has enough flush cycles,
- that only demineralised water is used for rinsing,
- that the air used to dry has been filtered (oil-free, low-microbiological contamination and particle-free) and,
- that the cleaning and disinfection unit is regularly serviced and tested.

The materials, concentrations, temperatures and treatment times, as well as rinsing requirements, specified by the cleaning and disinfection agent manufacturer must be adhered to at all times.

3.5.2.2 **Mechanical cleaning / disinfection (→ RECOMMENDED)**

Programme parameters used during certification (Miele G7835 CD cleaning and disinfection unit):

- 1 minute pre-cleaning (cold municipal water, drinking water quality <40°C) → Drain water → 3 minutes pre-cleaning (cold municipal water, drinking water quality <40°C) → Drain water
- 10 minutes cleaning at 55±5°C with 0.2% alkaline cleaning agent (0.2% Neodisher® MediClean) → Drain water
- 1 minute rinsing with demineralised water <40°C → Drain water → 2 minutes rinsing with demineralised water <40°C → Drain water
- Automatic disinfection > 5 minutes at 92±2°C with demineralised water.
- Automatic drying process 90±2°C of the cleaning and disinfecting unit for at least 30 minutes (± 60±5°C in the washing compartment).

(Re-)preparation process:

- Place the implements in the cleaning and disinfection unit. Make sure that the implements are not touching each other.
- Start the programme.
- When the programme ends, immediately remove the implements from the cleaning and disinfecting unit and ensure that they are dry enough before packaging.
- Inspect and package the implements as soon as possible after removing them from the unit.

3.5.2.3 **Manual subsequent drying**

If a manual subsequent drying is required, do so with a lint-free cloth and / or blow-out the lumens with sterile, oil-free pressurised air.

3.6 **Sterilisation – steam sterilisation – fractionated vacuum procedure**

Evidence of the fundamental suitability of the implement for an effective sterilisation was provided by an independent, accredited and ZLG-recognised (§15 (5) MGP) testing laboratory under use of an EHS3870 pre- and post-vacuum autoclave (Tuttnauer Europe B.V., Breda) and RB 51-3 and RB52-3 sterilised packaging (Steriking-foil). For this, the procedure described above was used.

3 vacuum cycles | 132°C / 270°F | ≥ 3 minutes stop time | Dry in the vacuum for at least 20 minutes³

3.6.1 Packaging

The sterilisation of the products must be done in suitable sterilisation packaging. The manufacturer's verification was done using doubled sterilisation packaging (hospital standard); that means that the single suitable sterilised foil packaging can also be used.

Flash sterilisation and the sterilisation of unpackaged implements is absolutely prohibited!

3.6.2 Sterilisation process – Fractionated vacuum procedure

Only the specified sterilisation procedures may be used for sterilisation.

Other sterilisation procedures are not permitted and their efficiency must be certified by the user / preparator themselves.

- **Fractionated vacuum procedure^{1,2} (with sufficient product drying³)**
- Steam steriliser conforming to DIN EN 13060/DIN EN 285 and ANSI AAMI ST 79 (FDA clearance in the USA)
- Validated in accordance with DIN EN ISO 17665 (valid IQ/OQ (commissioning) and product-specific performance evaluation (PQ))
- Maximum sterilisation temperature 134°C (273°F) including the tolerance margin in accordance with DIN EN ISO 17665
- Sterilisation time (time exposed to the sterilisation temperature)

4 Storage

After the sterilisation, the products must be stored dust-free and dry in the sterilisation packaging.

5 Repairs

¹ At least 3 vacuum stages

² The less effective gravitation procedure may only be used if the fractionated vacuum procedure is not available as an option, as it requires significantly longer sterilisation periods which must be determined and validated under the responsibility of the user, and specifically for their implements, devices, procedures and parameters.

³ The actual required product drying time depends directly on the parameters, which are solely the responsibility of the user (loading configuration and density, sterilisation status,...) and must therefore be determined by them. As a general rule, drying times should not be less than 20 minutes.

Josef Ganter GmbH provides a repair service for your implements (insofar as this is possible). If the implements are repaired by companies or persons which are not expressly authorised to do so by Josef Ganter, all guarantees shall become null and void.

Symbols used

Symbol	Description / Symbol for ...
	... Manufacturer
	... Article number
	... Production lot number – batch
	... is a medical device
	... Observe user instructions
	... CE conformity with the number of the notified body
	... Non-sterile
	... Warning

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Josef Ganter Feinmechanik GmbH